

## Quantitative RT-PCR for WNV and house-keeping gene ( $\beta$ -actin)

Primers :  
WNproC-F 10 (5'-CCTGTGTGAGCTGACAACTTAGT-3')  
WNproC-R 153 (5'-GCGTTTTAGCATATTGACAGCC-3')  
ACTBfrw966 : (5'-CAGCACAATGAAGATCAAGATCATC-3')  
ACTBrev1096 : (5'-CGGACTCATCGTACTCCTGCTT-3')

Resuspend primers with milliQ water to obtain a 100 $\mu$ M stock solution and store at -20°C (FMDV L3 laboratory).

### Probes :

Classical TaqMan probes :

WNproC-probe (5'-FAM-CCTGGTTTCTTAGACATCGAGATCT -TAMRA-3')  
ACTB1042-67 probe (5'-VIC-TCGCTGTCCACCTTCCAGCAGATGT-TAMRA-3')

Resuspend probes at 100 $\mu$ M with milliQ water (60 $\mu$ L for each probe). Prepare 5 $\mu$ L aliquotes and store at -20°C (FMDV L3 laboratory). Dilute probes at 10 $\mu$ M before use: add 45 $\mu$ L RNase free water for each 5  $\mu$ L aliquote.

### « Standard » :

• Viral RNA (strain IS98-ST1), stored at -80°C. Initial viral suspension titer: approximately 10<sup>6</sup>pfu/ $\mu$ L;

10-fold serial dilutions of viral RNA: 2 $\mu$ L previous dilution + 18 $\mu$ L RNase free water.

• extracted RNA from the horse « Deister »...

5-fold serial dilutions of equine RNA : 4 $\mu$ L previous dilution + 16 $\mu$ L RNase free water.

Kit : AgPath-ID™ One-Step RT-PCR Reagents (Applied Biosystem, 4387424 for 500 reactions), stored in -20°C freezer in Mix Preparation room in FMDV L3 laboratory.

PCR conditions: volumes per tube =

	ProC		ProC+ $\beta$ -actine
H <sub>2</sub> O DEPC	5.8 $\mu$ L	H <sub>2</sub> O DEPC	5.1 $\mu$ L
Tampon 2X	12,5 $\mu$ L	Tampon 2X	12,5 $\mu$ L
RT-PCR mix 25X	1 $\mu$ L	RT-PCR mix 25X	1 $\mu$ L
Primer ProC F 100 $\mu$ M	0.1 $\mu$ L	Amorce ProC F 100 $\mu$ M	0.1 $\mu$ L
Primer ProC R 100 $\mu$ M	0.1 $\mu$ L	Amorce ProC R 100 $\mu$ M	0.1 $\mu$ L
Probe ProC 10 $\mu$ M	0.5 $\mu$ L	Sonde ProC 10 $\mu$ M	0.5 $\mu$ L
		Amorce b-actin F 100 $\mu$ M	0.1 $\mu$ L
		Amorce b-actin R 100 $\mu$ M	0.1 $\mu$ L
		Sonde b-actin 10 $\mu$ M	0.5 $\mu$ L

Prepare a PCR mix, by calculating the volume required for the PCR reactions + 1 or 2 extra tubes. Distribute 20 $\mu$ L of PCR mix per PCR tube.

Add 5  $\mu$ L RNA or water per tube. Test dilutions 10<sup>-2</sup> to 10<sup>-6</sup> for the viral RNA (from 5x10<sup>4</sup> DECP<sub>50</sub> to 5 DECP<sub>50</sub>: one tube per dilution). Test dilutions 1/5 and 1/25 for equine RNA in the  $\beta$ -actin PCR. Add 1 NTC tube (5 $\mu$ L water) and 2 tubes per sample (5 $\mu$ L pure extracted RNA).

Enter the program on the AB7300 thermocycler (be careful to enter the right reaction volume = 25 $\mu$ L). Detector manager = FAM-TAMRA or VIC-TAMRA.

### Program :

Reverse transcription	45°C	10 min	
Denaturation	95°C	10 min	
PCR	95°C	15 s	} *40 cycles
Analysis mode = quantification	60°C	45s	